



A Publication
of Reliable Methods
for the Preparation
of Organic Compounds

Working with Hazardous Chemicals

The procedures in *Organic Syntheses* are intended for use only by persons with proper training in experimental organic chemistry. All hazardous materials should be handled using the standard procedures for work with chemicals described in references such as "Prudent Practices in the Laboratory" (The National Academies Press, Washington, D.C., 2011; the full text can be accessed free of charge at http://www.nap.edu/catalog.php?record_id=12654). All chemical waste should be disposed of in accordance with local regulations. For general guidelines for the management of chemical waste, see Chapter 8 of Prudent Practices.


In some articles in *Organic Syntheses*, chemical-specific hazards are highlighted in red "Caution Notes" within a procedure. It is important to recognize that the absence of a caution note does not imply that no significant hazards are associated with the chemicals involved in that procedure. Prior to performing a reaction, a thorough risk assessment should be carried out that includes a review of the potential hazards associated with each chemical and experimental operation on the scale that is planned for the procedure. Guidelines for carrying out a risk assessment and for analyzing the hazards associated with chemicals can be found in Chapter 4 of Prudent Practices.

The procedures described in *Organic Syntheses* are provided as published and are conducted at one's own risk. *Organic Syntheses, Inc.*, its Editors, and its Board of Directors do not warrant or guarantee the safety of individuals using these procedures and hereby disclaim any liability for any injuries or damages claimed to have resulted from or related in any way to the procedures herein.

These paragraphs were added in September 2014. The statements above do not supersede any specific hazard caution notes and safety instructions included in the procedure.

Organic Syntheses, Coll. Vol. 7, p.23 (1990); Vol. 61, p.98 (1983).

1,3-OXAZEPINES VIA PHOTOISOMERIZATION OF HETEROAROMATIC N-OXIDES: 3,1-BENZOXAZEPINE



Submitted by Angelo Albini, Gian Franco Bettinetti, and Giovanna Minoli¹.
Checked by Patrick MacManus and Robert M. Coates.

1. Procedure

Caution! 3,1-Benzoxazepine is a strong lacrimator and a moderate skin irritant. The preparation should be carried out in a well-ventilated hood. The apparatus should be shielded to avoid exposure to ultraviolet light.

Benzene has been identified as a carcinogen; OSHA has issued emergency standards on its use. All procedures involving benzene should be carried out in a well-ventilated hood, and glove protection is required.

Irradiation is carried out in a round-bottomed, cylindrical, Pyrex vessel (Note 1) equipped with a Pyrex immersion well (Note 2), nitrogen inlet, distillation sidearm, a small sidearm fitted with a rubber septum for removing aliquots, and a magnetic stirring bar. The flask is charged with 12 g (0.066 mol) (Note 3) of quinoline *N*-oxide dihydrate (Note 4) and 1.3 L of dry benzene (Note 5). The mixture is stirred and heated to boiling with a heating mantle as a slow stream of nitrogen (Note 6) is bubbled into the vessel. Benzene is distilled through the sidearm until the distillate is perfectly clear (Note 7). The lamp (Note 8) is placed in the immersion well, water is circulated through the cooling jacket, and the nitrogen flow is adjusted as necessary to maintain an outward flow of gas while the light yellow solution cools to room temperature (Note 9). The solution is stirred vigorously (Note 10) and irradiated for 2.5–3 hr, at which time the *N*-oxide is largely consumed (Note 11). The orange solution is transferred to a 1-L, round-bottomed flask, and the benzene is removed by rotary evaporation at room temperature. The red-orange oily residue, which contains some solid, is extracted with three 40-mL portions of dry cyclohexane, the combined cyclohexane extracts are evaporated under reduced pressure at room temperature, and the extraction operation is repeated on the oil thus obtained (Note 12). Evaporation of the combined cyclohexane extracts affords 6.1–6.3 g (63–65%) of crude 3,1-benzoxazepine (Note 13). Bulb-to-bulb distillation in a Kugelrohr apparatus at 0.2 mm with an oven temperature of 80°C affords 4.7–4.8 g (49–50%) of 3,1-benzoxazepine as a pale yellow oil, n_D^{24} 1.6074 (Note 14) and (Note 15).

2. Notes

1. The apparatus used by the checkers for irradiations at a 9-g scale was 33 cm high and 9 cm in diameter and had a 60/50 T joint at the top for the immersion well. The joints for the distillation and sampling sidearms were 24/40 and 14/20, respectively. The gas inlet was located about 6 cm from the bottom of the vessel to accommodate the use of a heating mantle. A disk of coarse, sintered glass was sealed into the gas inlet near its point of attachment. The capacity of the vessel with the immersion well in place was ca. 900 mL. The submitters used a similar but flat-bottomed apparatus of 1.2-L capacity. The flat-bottomed vessel facilitates vigorous stirring, but it does not fit as well into the heating mantle and may therefore be somewhat hazardous to use.

2. The checkers used a Vycor immersion well and a Pyrex filter sleeve. The immersion well, 450-W mercury lamp, and the requisite transformer are available from Hanovia Lamp Division, Canrad-Hanovia Inc., 100 Chestnut Street, Newark, NJ 07105.

3. The checkers carried out the irradiation on a 9-g scale in 750 mL of benzene after azeotropic

distillation.

4. [Quinoline *N*-oxide dihydrate](#) is supplied by Aldrich Chemical Company, Inc. and EGA Chemie KG, Steinheim/Albuch, Germany. The submitters prepared the compound by the procedure of Hayashi² with minor modifications. The water of hydration may be removed under reduced pressure in a drying pistol. However, since the anhydrous *N*-oxide is very hygroscopic, the submitters have found that it is more expedient to use the dihydrate and remove the water by azeotropic distillation in the irradiation vessel.

5. The checkers dried the [benzene](#) by distillation from [calcium hydride](#) immediately before use. The submitters report that [toluene](#) may be used instead of [benzene](#); however, since the product is not very thermally stable, they advise that the [toluene](#) should be evaporated without heating above 35–40°C during the isolation.

6. Other dry gases may be used. The submitters report that the reaction is not quenched by [oxygen](#).

7. A total of ca. 150–300 mL was collected. The distillation time may be reduced by insulating the vessel and sidearm with glass wool.

8. The submitters used a Helios Italquartz 500-W lamp which has emission characteristics similar to those of the Hanovia 450-W medium pressure [mercury](#) lamp used by the checkers.

9. The checkers noticed that a thin film of oil that was evidently [quinoline *N*-oxide](#) deposited on the surface of the immersion well and irradiation vessel during cooling.

10. Vigorous stirring is essential for optimum yields. The checkers obtained lower isolated yields (ca. 32–33%) in two runs in which a relatively slow stirring rate was employed. The low yields were probably caused in part by deposition of oil on the surface of the immersion well and the resulting interference with the transmission of UV light into the solution.

11. The submitters emphasize the importance of terminating the irradiation before all of the *N*-oxide is consumed. Overirradiation gives rise to a more complicated mixture of products from which the product can no longer be isolated by the simple extraction procedure described.

The progress of the irradiation was determined by the checkers by proton NMR analysis. At appropriate intervals 5-mL aliquots were removed, the solvent was evaporated, and [hexamethylbenzene](#) was added as an internal standard. The ratio of *N*-oxide, [benzoxazepine](#), and [hexamethylbenzene](#) was determined from integration of the resonances at δ 2.26 (s, 18H), 5.55 (d, 1H, $J = 6$), and 8.46 (d, 1H, $J = 6$), respectively, in [chloroform-*d*](#). After 2.5–3 hr of irradiation the amount of [benzoxazepine](#) present was ca. 60–68% of theoretical and ca. 10% of starting *N*-oxide remained.

The submitters followed the course of the irradiation by TLC analysis on silica gel with 5% (v/v) [methanol](#) in [chloroform](#) as developing solvent. Since some by-products have R_f values coincident with the *N*-oxide, this spot will not completely disappear and caution must be exercised to avoid overirradiation.

12. This extraction procedure separates most of the [carbostyryl](#), that is, [2\(1*H*\)-quinolinone](#), which is formed to the extent of ca. 20% in the irradiation. The submitters have isolated the [carbostyryl](#) by-product by crystallization of the extraction residue from 95% [ethanol](#) in runs carried out to high conversion. Alternatively, the [carbostyryl](#) may be isolated by chromatography of the crude product on silica gel with 5% [methanol-chloroform](#) as eluant. However, the [benzoxazepine](#) cannot be obtained by this method since it undergoes hydrolysis during the chromatography.

13. The purity of the crude product is about 90% according to NMR analysis, the remaining material being mostly unchanged *N*-oxide.

14. The spectral properties of the product are as follows: IR (liquid film) cm^{-1} : 1665 (C=N), 1630 (C=C), 1480 (sharp), 1440 (sharp), 1035 (strong), 765 (strong); ¹H NMR (CDCl₃) δ : 5.55 (d, 1 H, $J = 6$, CH=CH-O), 5.84 (d, 1 H, $J = 6$, CH=CH-O), 6.44 (s, 1 H, O-CH=N), 6.96 (m, 4 H, aromatic protons).


15. Samples of the product stored in tightly stoppered flasks in a freezer at –20°C for several weeks showed no sign of decomposition.

3. Discussion

The preparation of 3,1-benzoxazepines by photochemical isomerization of quinoline *N*-oxides constitutes a rather general entry into this class of seven-membered heterocycles. Since the structure of the photoisomer of [2-phenylquinoline *N*-oxide](#) was first recognized as [2-phenyl-3,1-benzoxazepine](#) by Buchardt et al.,³ the scope of this method for oxidative ring expansion of six-membered heterocyclic *N*-oxides to 1,3-oxazepines has been extensively explored.^{4 5} For example, irradiation of 2-cyano-, 2-phenyl-, and 2-methoxyquinoline *N*-oxides affords the corresponding 2-substituted 3,1-benzoxazepines

in 70–90% yield.^{6 7 8 9} However, isolation of the moisture-sensitive parent compound was only recently accomplished in the submitters' laboratories.¹⁰

Related 1,3-oxazepines have been obtained from irradiation of many other heterocyclic *N*-oxides including pyridine *N*-oxides, isoquinoline *N*-oxides, quinoxaline *N*-oxides, quinazoline *N*-oxides, phenanthridine *N*-oxides, benzophenazine *N*-oxides, and acridine *N*-oxides.^{4 5} However, the reported yields are variable and have generally been higher for phenyl and other aryl-substituted derivatives. This procedure is also satisfactory for 1,3-benzoxazepine (from isoquinoline *N*-oxide) and, in general, for oxazepines not substituent-stabilized.



A mechanism involving initial cyclization to an oxaziridine, [1,5] sigmatropic rearrangement to an imino epoxide, and electrocyclic ring opening was originally proposed for the photochemical isomerization.^{4 5} However, since later attempts to detect intermediates by flash photolysis were unsuccessful,^{11 12} ground-state oxaziridines, if formed at all, must have exceedingly short lifetimes. The benzoxazepines undergo facile hydrolysis to *o*-(*N*-acylamino)phenylacetaldehydes, which frequently exist as the cyclic carbinol amide tautomers. If water is present during the irradiation from use of the *N*-oxide hydrate or moist solvent, the hydrolysis products may be isolated instead of the benzoxazepine. Dehydration of the carbinol amide to *N*-acyl indoles may also occur during irradiation and/or purification of the products. The formation of carbostyrils is sometimes an important competing reaction in the irradiation of quinoline *N*-oxides, and this by-product is, in fact, formed to the extent of ca. 20% in the present procedure. The use of polar protic solvents such as water or alcohols favors carbostyryl formation in contrast to aprotic solvents such as benzene or acetone in which the pathway leading to benzoxazepines usually predominates.

References and Notes

1. Istituto di Chimica Organica dell'Università, v. le, Taramelli 10, 27100 Pavia, Italy.
 2. Hayashi, E., private communication to Ochiai, E. in "Aromatic Amine *N*-Oxides," Elsevier Publishing Co.: Amsterdam, 1967, p. 24.
 3. Buchardt, O.; Jensen, B.; Larsen, I. K. *Acta Chem. Scand.* **1967**, *21*, 1841–1854.
 4. Spence, G. G.; Taylor, E. C.; Buchardt, O. *Chem. Rev.* **1970**, *70*, 231–265;
 5. Bellamy, F.; Streith, J. *Heterocycles* **1976**, *4*, 1391–1447.
 6. Kaneko, C.; Yamada, S. *Chem. Pharm. Bull.* **1966**, *14*, 555–557;
 7. Kaneko, C.; Yamada, S.; Ishikawa, M. *Tetrahedron Lett.* **1966**, 2145–2150;
 8. Buchardt, O.; Kumler, P. L.; Lohse, C. *Acta Chem. Scand.* **1969**, *23*, 1155–1167
 9. Albini, A.; Fasani, E.; Dacrema, L. M. *J. Chem. Soc. Perkin Trans. 1* **1980**, 2738–2742.
 10. Albini, A.; Bettinetti, G. F.; Minoli, G. *Tetrahedron Lett.* **1979**, 3761–3764.
 11. Lohse, C. *J. Chem. Soc., Perkin Trans. II* **1972**, 229–233;
 12. Tomer, K. B.; Harrit, N.; Rosenthal, I.; Buchardt, O.; Kumler, P. L.; Creed, D. *J. Am Chem. Soc.* **1973**, *95*, 7402–7406.
-

Appendix

**Chemical Abstracts Nomenclature (Collective Index Number);
(Registry Number)**

silica gel

1,3-OXAZEPINES

3,1-benzoxazepines

quinoline N-oxides

pyridine N-oxides

isoquinoline N-oxides

quinoxaline N-oxides

quinazoline N-oxides

phenanthridine N-oxides

benzophenazine N-oxides

acridine N-oxides

carbostyrils

benzoxazepines

[ethanol \(64-17-5\)](#)

[Benzene \(71-43-2\)](#)

[methanol \(67-56-1\)](#)

[chloroform \(67-66-3\)](#)

[oxygen \(7782-44-7\)](#)

[nitrogen \(7727-37-9\)](#)

[mercury \(7439-97-6\)](#)

[cyclohexane \(110-82-7\)](#)

[acetone \(67-64-1\)](#)

[toluene \(108-88-3\)](#)

[Hexamethylbenzene \(87-85-4\)](#)

carbostyryl (59-31-4)

calcium hydride (7789-78-8)

quinoline N-oxide (1613-37-2)

isoquinoline N-oxide (1532-72-5)

chloroform-d (865-49-6)

3,1-Benzoxazepine (15123-59-8)

benzoxazepine

2-phenyl-3,1-benzoxazepine

1,3-benzoxazepine

quinoline N-oxide dihydrate

2(1H)-quinolinone

2-phenylquinoline N-oxide