



A Publication
of Reliable Methods
for the Preparation
of Organic Compounds

Working with Hazardous Chemicals

The procedures in *Organic Syntheses* are intended for use only by persons with proper training in experimental organic chemistry. All hazardous materials should be handled using the standard procedures for work with chemicals described in references such as "Prudent Practices in the Laboratory" (The National Academies Press, Washington, D.C., 2011; the full text can be accessed free of charge at http://www.nap.edu/catalog.php?record_id=12654). All chemical waste should be disposed of in accordance with local regulations. For general guidelines for the management of chemical waste, see Chapter 8 of Prudent Practices.


In some articles in *Organic Syntheses*, chemical-specific hazards are highlighted in red "Caution Notes" within a procedure. It is important to recognize that the absence of a caution note does not imply that no significant hazards are associated with the chemicals involved in that procedure. Prior to performing a reaction, a thorough risk assessment should be carried out that includes a review of the potential hazards associated with each chemical and experimental operation on the scale that is planned for the procedure. Guidelines for carrying out a risk assessment and for analyzing the hazards associated with chemicals can be found in Chapter 4 of Prudent Practices.

The procedures described in *Organic Syntheses* are provided as published and are conducted at one's own risk. *Organic Syntheses, Inc.*, its Editors, and its Board of Directors do not warrant or guarantee the safety of individuals using these procedures and hereby disclaim any liability for any injuries or damages claimed to have resulted from or related in any way to the procedures herein.

September 2014: The paragraphs above replace the section "Handling and Disposal of Hazardous Chemicals" in the originally published version of this article. The statements above do not supersede any specific hazard caution notes and safety instructions included in the procedure.

Copyright © 2007 Organic Syntheses, Inc. All Rights Reserved

ELECTROPHILIC CYCLIZATION WITH *N*-IODOSUCCINIMIDE: PREPARATION OF 5-(4-BROMOPHENYL)-3-iodo-2-(4-METHYL- PHENYL)FURAN



Submitted by Adam Sniady, Marco S. Morreale, and Roman Dembinski.¹

Checked by Cody Timmons and Peter Wipf.²

1. Procedure

A. 1-(4-Bromophenyl)-4-(4-methylphenyl)but-3-yn-1-ol (1). A 500-mL, three-necked, round-bottomed flask (Note 1), equipped with a magnetic stirring bar and fitted with a nitrogen inlet adapter and two rubber septa, is charged with THF (40 mL) (Note 2) and 4-ethynyltoluene (4.08 mL, 32.2 mmol) (Note 3) via syringes. The reaction vessel is placed in a dry ice-acetone bath. After 10 min, lithium diisopropylamide solution (LDA) (16.1 mL, 32.2 mmol) (Note 4) is added within 5 min via syringe to the stirred reaction mixture. The mixture is allowed to stir for 10 min after which time the bath is removed and the reaction vessel is allowed to warm to room temperature. Dimethyl sulfoxide (DMSO) (75 mL) (Note 5) is added over 5 min via syringe to the reaction vessel, followed by a solution of 4-bromostyrene oxide (4.94 g, 24.8 mmol) (Note 6) in 5 mL of THF over 1 min via syringe. The reaction mixture acquires a dark red wine color and is allowed to stir at room temperature for 2 h (Note 7). Then, the reaction flask is placed in an ice-bath and quenched with saturated aqueous ammonium

chloride (100 mL). The mixture is transferred to a 500-mL separatory funnel using saturated aqueous ammonium chloride (50 mL) as a rinse. The aqueous phase is extracted with ether (3×150 mL). The combined ether extracts are transferred to a 1-L separatory funnel and washed with saturated aqueous sodium chloride (2×200 mL). The organic phase is dried with anhydrous sodium sulfate (10 g) and filtered through a fritted-glass funnel of medium porosity. The solvent is removed via rotary evaporation (23 °C, 25 mmHg) followed by high vacuum (23 °C, 0.6 mmHg, 3 h). The crude product is purified by recrystallization and chromatography of the mother liquor on SiO₂ (Note 8) to yield **1** as a white crystalline solid (6.01 g, 19.1 mmol, 77%) (Note 9).

B. 1-(4-Bromophenyl)-4-(4-methylphenyl)but-3-yn-1-one (2). A 500-mL, single-necked, round-bottomed flask equipped with a rubber septum and a magnetic stirring bar is charged with **1** (6.01 g, 19.1 mmol) and acetone (200 mL) (Note 10). Jones reagent (8.0 mL, 3.0 M solution, 24 mmol) (Note 11) is added by syringe over 5 min at 0 °C (ice-water bath) to the stirred mixture (Notes 12, 13). The bath is removed and the reaction mixture is allowed to stir for an additional 15 min (Note 7), then immediately ether (130 mL) and isopropanol (4 mL) are added. The blue-green solid is filtered off through a medium porosity fritted-glass funnel and the filtrate is transferred to a 1-L separatory funnel filled with 100 mL of ice water. Three extractions with ether (3×100 mL) are carried out. The combined organic phases are washed twice with saturated sodium bicarbonate (2×150 mL), dried over sodium sulfate (10 g), filtered, and concentrated at reduced pressure on a rotary evaporator (23 °C, 25 mmHg). The resulting light yellow solid is transferred to a 50-mL, single neck round-bottomed flask and dried for 6 h under vacuum (23 °C, 0.6 mmHg) to give **2** (5.56 g, 17.8 mmol, 93%), which is used in the next step without further purification (Note 14).

C. 5-(4-Bromophenyl)-3-iodo-2-(4-methylphenyl)furan (3). A 500-mL, three-necked, round-bottomed flask equipped with a magnetic stirring bar and a rubber septum is charged with **2** (5.10 g, 16.3 mmol). An additional rubber septum with a nitrogen inlet is attached to the remaining neck, and acetone (150 mL) (Note 10) is added via syringe. In a separate flask *N*-iodosuccinimide (NIS) (3.85 g, 17.1 mmol) (Note 15) is dissolved in acetone (60 mL). The solution of NIS is then injected via syringe over 5 min into the reaction mixture that is stirred at room temperature under a nitrogen atmosphere for 2 h. TLC analysis (Note 7) indicates a complete

conversion. The solvent is removed by rotary evaporation (23 °C, 25 mmHg) and the residue is extracted with ether (150 mL) (Note 16). The solid (succinimide) is filtered off through a medium porosity fritted-glass funnel, and additionally washed with ether (60 mL) (Note 16). The ethereal extracts are combined and transferred to a 500-mL separatory funnel containing 50 mL of a 2.0 M aqueous sodium bisulfite solution. The phases are separated and the aqueous phase is further extracted with ether (3 × 50 mL). The combined organic extracts are washed with brine (1 × 100 mL) and dried over Na₂SO₄ (10 g). The mixture is filtered through a medium porosity fritted glass funnel and the solvent is removed by rotary evaporation (23 °C, 25 mmHg). The resulting material is charged on a column (3 × 40 cm) of silica gel (35 g) (Note 17) and eluted with hexanes/ethyl acetate (8:1) (approximately 600 mL) (Note 18). The product is concentrated by rotary evaporation (23 °C, 25 mmHg) and then dried at 0.6 mmHg (23 °C, 6 h) to yield **3** as a colorless solid (5.80 g, 13.2 mmol, 81%) (Note 19).

2. Notes

1. The apparatus was oven-dried at 135 °C for 2 h and maintained under an atmosphere of dry nitrogen during the course of the reaction. The flasks had neck sizes of 24/40.

2. THF was dried over sodium/benzophenone and freshly distilled under an atmosphere of dry nitrogen.

3. 4-Ethynyltoluene (98%) was purchased from GFS Chemicals and used as received. The checkers used chemicals from Alfa Aesar without further purification.

4. Lithium diisopropylamide (2.0 M solution in THF/*n*-heptane) was purchased from Acros Organics and used as received.

5. Methylsulfoxide, anhydrous, DriSolv was purchased from EMD Chemicals, Inc. and used as received. The checkers used 99.7% anhydrous DMSO from Acros Organics.

6. 4-Bromostyrene oxide (96%) was purchased from Aldrich Chemical Company, Inc. and used as received.

7. The progress of the reaction was monitored by TLC analysis on silica gel eluted with hexanes/ethyl acetate (8:2) and visualized with phosphomolybdic acid. The *R_f* values were: 0.53 (blue) for 4-bromostyrene oxide, 0.32 (dark blue) for alcohol **1**, 0.59 (blue) for ketone **2**, and 0.74 (orange-brown) for furan **3**.

8. Ethyl acetate (15 mL) was added to the flask containing crude **1** (10.23 g). The flask was heated gently with a heat-gun under hand-stirring until the all of the crude material was dissolved. Then, the solution was layered with hexanes (30 mL) and placed in the refrigerator (−6 °C) for 12 h to induce crystallization. The flask was placed in an ice-bath and the mother liquor was removed via syringe. The crystalline residue was rinsed with cold hexanes (2 × 10 mL) (placed in an ice bath for at least 45 min prior to use) and the rinses were removed via syringe; then the crystalline material was dried (23 °C, 0.6 mmHg, 3 h) to afford alcohol **1** (4.45 g, 14.1 mmol, 57%). The rinses were combined with the mother liquor and the solvent was removed via rotary evaporation (23 °C, 25 mmHg) followed by high vacuum (23 °C, 0.6 mmHg, 3 h) to yield 4.12 g of crude material. This material was dissolved in ethyl acetate (10 mL) with hand-stirring and the help of a heat-gun. The solution was layered with hexanes (20 mL) and placed in the refrigerator for 48 h. Then the crystallization flask was placed into an ice bath and the mother liquor was removed via syringe. The crystalline material was rinsed with cold hexanes (10 mL) and the rinse was removed via syringe; then the crystalline material was dried (23 °C, 0.6 mmHg, 3 h) to afford an additional 6% of alcohol **1** (0.46 g, 1.5 mmol). The rinses and mother liquor were combined and concentrated by rotary evaporation. The resulting material was charged on a column (3 × 40 cm) of silica gel (35 g) and eluted with hexanes/ethyl acetate (6:1) (approximately 400 mL) to provide an additional 14% of alcohol **1** (1.10 g, 3.49 mmol). The total amount of alcohol **1** obtained in this reaction was 6.01 g (19.1 mmol, 77%).

9. mp 133–135 °C (ethyl acetate/hexanes). Calcd for C₁₇H₁₅BrO: C, 64.78; H, 4.80. Found: C, 64.59; H, 4.85. IR (cm^{−1}, NaCl disk) 3411 (m), 1509 (m), 1046 (m), 817 (s). UV-Vis (ε, M^{−1}cm^{−1}; ether; 3.5 × 10^{−5} M) 244 (29000), 254 (24000). MS (rel intensity) 185 (88), 130 (100), 115 (74), 77 (96), 63 (34). ¹H NMR (300 MHz, CDCl₃) δ: 2.36 (s, 3 H), 2.61 (d, 1 H, *J* = 3.0 Hz), 2.80–2.84 (m, 2 H), 4.90 (td, 1 H, *J* = 6.3, 1.5 Hz), 7.12 (d, 2 H, *J* = 7.8 Hz), 7.28–7.33 (m, 4 H), 7.47 (d, 2 H, *J* = 8.4 Hz); ¹³C NMR δ: 21.4, 30.6, 71.9, 83.6, 84.5, 119.9, 121.6, 127.5, 129.0, 131.4, 131.5, 138.2, 141.6.

10. High Purity Solvent grade purchased from Burdick & Jackson was used as received. The checkers used OmniSolv High Purity Solvent grade acetone from EMD Chemicals. For procedure B reagent grade acetone can be used.

11. Prepared according to ref. 11.

12. Precipitation of blue-green chromium(III) salts was observed.

13. Reaction can be carried out in an open flask.

14. The light yellow solid had a purity of approximately 95% as determined by NMR. A portion of the material (0.46 g) was retained for recrystallization as described in Note 8 to afford colorless crystals for analytical data: mp 121–123 °C (ethyl acetate /hexanes). Calcd for C₁₇H₁₃BrO: C, 65.19; H, 4.18. Found: C, 65.33; H, 4.14. IR (cm⁻¹, NaCl disk) 1677 (s), 809 (s). UV-Vis (ε, M⁻¹cm⁻¹; ether; 4.5 × 10⁻⁵ M) 254 (35000). MS (rel intensity) 312 (M⁺, 53), 183 (100), 155 (60), 129 (74). ¹H NMR (CDCl₃) δ: 2.33 (s, 3 H), 4.02 (s, 2 H), 7.09 (d, 2 H, *J* = 7.8 Hz), 7.31 (d, 2 H, *J* = 8.1 Hz), 7.61 (d, 2 H, *J* = 8.4 Hz), 7.91 (d, 2 H, *J* = 8.4 Hz); ¹³C NMR δ: 21.3, 31.5, 80.9, 85.7, 119.7, 128.7, 128.9, 130.1, 131.4, 131.9, 134.1, 138.3, 192.1.

15. NIS was obtained from Simafex and used as received. The checkers used NIS from Alfa Aesar.

16. Ether (unstabilized) was purchased from EMD Chemicals, Inc. and used as received. The checkers used diethyl ether that was freshly distilled from sodium/benzophenone. Extraction with ether stabilized with ethanol led to “leaching” of the succinimide.

17. Sorbent Technologies 60 Å, 230-400 mesh.

18. Flash column chromatography was performed in order to remove residues of succinimide. The solid was dissolved in a minimum amount of hexanes/ethyl acetate (8:1) and loaded onto the column. Pressured air was applied and one fraction was collected.

19. mp 98–100 °C (ethyl acetate /hexanes). Calcd for C₁₇H₁₂BrIO: C, 46.50; H, 2.75. Found: C, 46.72; H, 2.74. IR (cm⁻¹, NaCl disk) 1492 (s), 1475 (s), 1400 (w), 946 (s), 816 (s). UV-Vis (ε, M⁻¹cm⁻¹; ether, 2.5 × 10⁻⁵ M) 237 (23000), 253 (22000), 330 (35000). MS (rel intensity) 438 (M⁺, 100), 283 (M – C₆H₄Br, 70). ¹H NMR (CDCl₃) δ: 2.42 (s, 3 H), 6.83 (s, 1 H), 7.28 (d, 2 H, *J* = 8.1 Hz), 7.48–7.59 (m, 4 H), 7.94 (d, 2 H, *J* = 8.1 Hz); ¹³C NMR δ: 21.4, 62.1, 116.1, 121.7, 125.2, 126.1, 127.2, 128.5, 129.1, 131.9, 138.4, 151.7, 152.6.

Safety and Waste Disposal Information

All hazardous materials should be handled and disposed of in accordance with “Prudent Practices in the Laboratory”; National Academy Press; Washington, DC, 1995.

3. Discussion

The furan ring is a key component of many pharmaceutical agents, fragrances, natural products, and organic materials. Halofurans allow for further functionalization of the ring, usually via transition metal-catalyzed coupling reactions; they serve as starting materials for combinatorial chemistry.³ The title molecule includes differently halogenated furan and phenyl rings that can be further functionalized.


The method described herein applies electrophilic cyclization for the synthesis of 3- (*beta*) iodofurans, with different, but electronically similar substituents in the 2- and 5- (*alpha*) positions. The method is of general use; other 3-halo (iodo- and bromo-) furans can also be synthesized in a regioselective manner (Table 1),⁴ using other electrophilic reagents (ICl, NBS, and Br₂). This approach provides a facile, non-acidic, non-basic, non-transition metal based protocol, in which an easy to handle reagent⁵ facilitates cyclization and serves simultaneously as a halogen donor. The reaction conditions are mild (room temperature in acetone). The yield is high, the reaction time is relatively short (0.5–2 h), and the isolation procedure is simple.

The best current alternatives for the synthesis of 2,5-unsymmetrically substituted 3-iodofurans include: a) Hydriodic acid cyclization of 4-hydroxybut-2-yn-1-ones.^{3,6} This procedure features a strong acid. b) Iodocyclization (3 equiv of iodine) of alk-3-yn-1,2-diols that requires the use of a base.⁷ c) Sonogashira coupling-addition-cyclocondensation one-pot procedure that starts from acyl chlorides and tetrahydropyranyl ethers of propargyl alcohols.⁸

Alkynol **1** was prepared by a known method.⁹ It was observed by ¹H NMR that the ring opening of the 2-(4-bromophenyl)oxirane proceeded in a completely regioselective manner. BF₃·Et₂O was not used as an additive, as referenced; although such a procedure may be carried out in THF.¹⁰ THF and DMF were found to be inefficient solvents for this protocol. Alkynone **2** was prepared by oxidation of the alkynol **1** with CrO₃ (Jones reagent).¹¹ Among other oxidation methods screened, Dess-Martin periodinane¹² was comparable. Jones oxidation for relevant compounds has been reported with low,¹³ reasonable,^{14,15} and high yields.¹⁶ No experimental procedure, however, was provided for the high yielding example. A catalytic oxidation protocol has also been described.¹⁷ Alk-3-yn-1-ones are reportedly prone to isomerization to the corresponding allenes;^{17,18} this does not pertain to the

procedure for **2**, presumably due to stabilization of the triple bond by conjugation to the phenyl ring. Alkynone **2** was found to be stable as a solid (stored in a freezer); cycloisomerization was observed when **2** was left for a prolonged time in solution. Cycloisomerization of homologous alk-3-yn-1-ones to the non-halogenated furans has been reported to proceed efficiently with the use of metal catalysts such as Hg,¹⁶ Au,¹⁹ Pd,^{15,20} or Zn²¹ derivatives.

Table 1. Preparation of other halofurans via electrophilic cyclization.⁴



<i>Entry</i>	<i>R</i>	<i>R'</i>	<i>E-A, solvent</i>	<i>Yield (%)</i>
1	C ₆ H ₅	<i>p</i> -CH ₃ C ₆ H ₄	NBS, acetone	89
2			NIS, acetone	84
3			ICl, CH ₂ Cl ₂	82
4	<i>p</i> -BrC ₆ H ₄	<i>p</i> -CH ₃ C ₆ H ₄	NBS, acetone	89
5			NIS, acetone	86
6	<i>p</i> -ClC ₆ H ₄	<i>p</i> -CH ₃ C ₆ H ₄	NBS, acetone	83
7			NIS, acetone	94
8	<i>p</i> -BrC ₆ H ₄	<i>p-t</i> -BuC ₆ H ₄	NBS, acetone	87
9			NIS, acetone	81
10	<i>p</i> -ClC ₆ H ₄	<i>p-t</i> -BuC ₆ H ₄	NBS, acetone	86
11			NIS, acetone	88

1. Department of Chemistry, Oakland University, Rochester, MI 48309-4477, USA; dembinsk@oakland.edu
2. Department of Chemistry, University of Pittsburgh, Pittsburgh, PA 15260, USA; pwipf@pitt.edu
3. Sauers, R. R.; Van Arnum, S. D. *J. Comb. Chem.* **2004**, *6*, 350.

4. For reaction leading to 3-chlorofurans see: Sniady, A.; Morreale, M. S.; Wheeler, K. A.; Dembinski, R. *Eur. J. Org. Chem.* **2008**, 3449.
5. Sniady, A. *Synlett* **2006**, 960.
6. Obrecht, D. *Helv. Chim. Acta* **1989**, 72, 447.
7. a) El-Taeb, G. M. M.; Evans, A. B.; Jones, S.; Knight, D. W. *Tetrahedron Lett.* **2001**, 42, 5945. b) Bew, S. P.; El-Taeb, G. M. M.; Jones, S.; Knight, D. W.; Tan, W.-F. *Eur. J. Org. Chem.* **2007**, 5759.
8. a) Karpov, A. S.; Merkul, E.; Oeser, T.; Müller, T. J. *J. Chem. Commun.* **2005**, 2581. b) Karpov, A. S.; Merkul, E.; Oeser, T.; Müller, T. J. *Eur. J. Org. Chem.* **2006**, 2991.
9. Brandsma, L. *Preparative Acetylenic Chemistry*, 2nd Ed.; Studies in Organic Chemistry 34; Elsevier: Amsterdam, 1988; p 67.; a relevant procedure that uses ethylenediamine complex of lithium acetylide is also available: McDonald, F. E.; White, B. H. *Org. Synth. Coll. Vol.* 10, 621 (2004); Vol. 79, 27 (2002).
10. Yamaguchi, M.; Hirao, I. *Tetrahedron Lett.* **1983**, 24, 391.
11. Eisenbraun, E. *J. Org. Synth. Coll. Vol.* 5, 310 (1973); Vol. 45, 28 (1965).
12. a) Boeckman, Jr., R. K.; Shao, P. ; Mullins, J. J. *Org. Synth. Coll. Vol.* 10, 696 (2004); Vol. 77, 141 (2000). b) Ireland, R. E.; Liu, L. *J. Org. Chem.* **1993**, 58, 2899.
13. Bol'shedvorskaya, R. L.; Pavlova, G. A.; Gavrilov, L. D.; Alekseeva, N. V.; Vereshchagin, L. I. *Zh. Org. Khim.* **1972**, 8, 1879; *J. Org. Chem. USSR* **1972**, 8, 1927.
14. Padwa, A.; Crumrine, D.; Hartman, R.; Layton, R. *J. Am. Chem. Soc.* **1967**, 89, 4435.
15. Fukuda, Y.; Shiragami, H.; Utimoto, K.; Nozaki, H. *J. Org. Chem.* **1991**, 56, 5816.
16. Ménard, D.; Vidal, A.; Barthomeuf, C.; Lebreton, J.; Gosselin, P. *Synlett* **2006**, 57.
17. Schmieder-van de Vondervoort, L.; Bouttemy, S.; Padrón, J. M.; Le Bras, J.; Muzart, J.; Alsters, P. L. *Synlett* **2002**, 243.
18. Krause, N.; Hashmi, A. S. K. *Modern Allene Chemistry*, Wiley: Weinheim, 2004; p 3-50.
19. Hashmi, A. S. K.; Schwarz, L.; Choi, J.-H.; Frost, T. M. *Angew. Chem., Int. Ed.* **2000**, 39, 2285.
20. a) Utimoto, K. *Pure Appl. Chem.* **1983**, 55, 1845. b) Sheng, H.; Lin, S.; Huang, Y. *Synthesis* **1987**, 1022.

21. a) Sniady, A.; Durham, A.; Morreale, M. S.; Wheeler, K. A.; Dembinski, R. *Org. Lett.* **2007**, 9, 1175. b) Sniady, A.; Durham, A.; Morreale, M. S.; Marcinek, A.; Szafert, S.; Lis, T.; Brzezinska, K. R.; Iwasaki, T.; Ohshima, T.; Mashima, K.; Dembinski, R. *J. Org. Chem.* **2008**, 73, 5881.

Appendix

Chemical Abstracts Nomenclature; (Registry Number)

5-(4-Bromophenyl)-3-iodo-2-(4-methylphenyl)furan; (852211-52-0)

2-(4-Bromophenyl)oxirane: (4-Bromostyrene oxide); (32017-76-8)

4-Ethynyltoluene; (9766-97-2)

Lithium diisopropylamide; (4111-54-0)

1-(4-Bromophenyl)-4-(4-methylphenyl)but-3-yn-1-ol; (852211-40-6)

Chromium trioxide; (1333-82-0)

1-(4-Bromophenyl)-4-(4-methylphenyl)but-3-yn-1-one; (852211-45-1)

N-Iodosuccinimide (NIS); (516-12-1)



Roman Dembinski was educated at the Polish Academy of Sciences, Centre of Molecular and Macromolecular Studies, Łódź, Poland (Ph.D. with A. Skowrońska). Before, he completed his undergraduate studies at the Łódź University of Technology. He then joined the research group of J. A. Gladysz at the University of Utah as a research associate. Since 1999 he has held a faculty appointment at the Oakland University, Rochester. His scientific interests include synthetic organic, fluororous, medicinal, and organometallic chemistry.



Adam Sniady obtained his M.S. in chemistry from the Łódź University of Technology, Poland, in 2004. Earlier, he received M.S. in Biology from the University of Łódź. He is currently pursuing a Ph.D. with Roman Dembinski at the Oakland University, Rochester, Michigan, USA. His research is focused on developing new synthetic methodologies, including electrophilic cyclizations and synthesis of modified nucleosides.




Marco Morreale obtained his B.A. in chemistry from Wayne State University (Detroit, Michigan) in 1999. During his undergraduate years, his introduction to chemical research began in the laboratory of Norman Lebel. He entered the graduate program in chemistry at Oakland University (Rochester, Michigan) in the fall of 2004. Under the direction of Roman Dembinski, he is currently investigating cyclization strategies for the synthesis of substituted furans.



Cody Timmons received his B.S. (2000) and M.S. (2002) degrees in Chemistry from West Texas A&M University in Canyon, Texas. His M.S. research was performed under the direction of Professor Darren Williams and dealt with the kinetics of the degradation of explosives. He then moved to Texas Tech University in the fall of 2002 and joined the group of Professor Guigen Li, where he studied novel olefin addition and carbon-carbon bond formation methodologies. He completed his Ph.D. in 2006 and is currently a postdoctoral research associate in the lab of Professor Peter Wipf at the University of Pittsburgh, where he is pursuing the total synthesis of ergot and clavine alkaloids.

CT060, CDCl3, 75MHz




Current Data Parameters
NAME CT6013C
EXPNO 3
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060722
Time 15.53
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zgpg
TD 65536
SOLVENT CDCl3
NS 32
DS 4
SWH 18115.941 Hz
FIDRES 0.276427 Hz
AQ 1.8088436 sec
RG 4597.6
DW 27.600 usec
DE 6.00 usec
TE 300.0 K
D1 10.00000000 sec
d11 0.03000000 sec
DELTA 9.89999962 sec
TD0 1

===== CHANNEL f1 =====
NUC1 13C
P1 9.00 usec
PL1 1.80 dB
SFO1 75.5381641 MHz

===== CHANNEL f2 =====
CPDPRG2 waltz16
NUC2 1H
PCPD2 92.00 usec
PL2 0.00 dB
PL12 20.00 dB
PL13 20.00 dB
SFO2 300.3812015 MHz

F2 - Processing parameters
SI 65536
SF 75.5306166 MHz
WDW EM
SSB 0
LB 1.00 Hz
GB 0
PC 1.40



CT104, 75 MHz, ¹³C, CDC13




Current Data Parameters
NAME CT104C13
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060901
Time 8.19
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zgpg
TD 65536
SOLVENT CDC13
NS 64
DS 4
SWH 18115.941 Hz
FIDRES 0.276427 Hz
AQ 1.8088436 sec
RG 10321.3
DW 27.600 usec
DE 6.00 usec
TE 300.0 K
D1 6.00000000 sec
d11 0.03000000 sec
DELTA 5.90000010 sec
TD0 1

===== CHANNEL f1 =====
NUC1 ¹³C
P1 9.00 usec
PL1 1.80 dB
SFO1 75.5381641 MHz

===== CHANNEL f2 =====
CPDPRG2 waltz16
NUC2 ¹H
PCPD2 92.00 usec
PL2 0.00 dB
PL12 20.00 dB
PL13 20.00 dB
SFO2 300.3812015 MHz

F2 - Processing parameters
SI 65536
SF 75.5306208 MHz
WDW EM
SSB 0
LB 1.00 Hz
GB 0
PC 1.40



furan3, 13C, CDC13, 75MHz




Current Data Parameters
NAME furan3C
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060906
Time_ 18.50
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zgpg
TD 65536
SOLVENT CDC13
NS 192
DS 4
SWH 18115.941 Hz
FIDRES 0.276427 Hz
AQ 1.8088436 sec
RG 8192
DW 27.600 usec
DE 6.00 usec
TE 300.0 K
D1 5.00000000 sec
d11 0.03000000 sec
DELTA 4.90000010 sec
TD0 1

===== CHANNEL f1 =====
NUC1 13C
P1 9.00 usec
PL1 1.80 dB
SFO1 75.5381641 MHz

===== CHANNEL f2 =====
CPDPRG2 waltz16
NUC2 1H
PCPD2 92.00 usec
PL2 0.00 dB
PL12 20.00 dB
PL13 20.00 dB
SFO2 300.3812015 MHz



F2 - Processing parameters
SI 65536
SF 75.5306147 MHz
WDW EM
SSB 0
LB 1.00 Hz
GB 0
PC 1.40



furan3, 1H, CDC13, 300MHz

7.956
7.929
7.563
7.537
7.289
7.270
7.262
6.832

2.419




Current Data Parameters
NAME furan3H
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060906
Time_ 18.22
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDC13
NS 8
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 114
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799985 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



furan3, 1H, CDC13, 300MHz




Current Data Parameters
NAME furan3H
EXPNO 1
PROCNO 1


F2 - Acquisition Parameters
Date_ 20060906
Time_ 18.22
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDC13
NS 8
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 114
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799985 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT060, CDC13, 300MHz




Current Data Parameters
NAME CT60C13
EXPNO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060722
Time_ 15.35
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDC13
NS 4
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 71.8
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799995 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT060, CDCl3, 300MHz




Current Data Parameters
NAME CT60C13
EXPNO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060722
Time 15.35
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDCl3
NS 4
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 71.8
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799995 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT060, CDC13, 300MHz




Current Data Parameters
NAME CT60C13
EXPNO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060722
Time_ 15.35
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDC13
NS 4
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 71.8
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799995 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

4.923
4.909
4.904
4.891
4.883



CT060, CDC13, 300MHz




Current Data Parameters
NAME CT60C13
EXPNO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060722
Time 15.35
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDC13
NS 4
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 71.8
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799995 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT104, 300 MHz, 1H, CDCl3




Current Data Parameters
NAME CT104H
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060901
Time_ 8.03
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDCl3
NS 8
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 35.9
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799993 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT104, 300 MHz, 1H, CDCl3




Current Data Parameters
NAME CT104H
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060901
Time_ 8.03
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDCl3
NS 8
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 35.9
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799993 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



CT104, 300 MHz, 1H, CDCl3



Current Data Parameters
NAME CT104H
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20060901
Time_ 8.03
INSTRUM spect
PROBHD 5 mm DUL 1H-13
PULPROG zg
TD 65536
SOLVENT CDCl3
NS 8
DS 2
SWH 6218.905 Hz
FIDRES 0.094893 Hz
AQ 5.2691445 sec
RG 35.9
DW 80.400 usec
DE 6.00 usec
TE 300.0 K
D1 2.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 9.00 usec
PL1 1.00 dB
SFO1 300.3818550 MHz

F2 - Processing parameters
SI 32768
SF 300.3799993 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

