**Date:**

4 Nov 2009

**Document's Title:**

fid

**Spectrum Title:**

User/Group Lars Troendlin/Pfaltz

LT445 nach Auf.

\_Proton CDCl3 uicon troend 14

**Frequency (MHz):**

(f1) 400.132

**Original Points Count:**

(f1) 24576

**Actual Points Count:**

(f1) 65536

**Acquisition Time (sec):**

(f1) 2.9688

**Spectral Width (ppm):**

(f1) 20.689

**Pulse Program:**

ZG30

**Temperature:**


300

**Number of Scans:**

16

**Acq. Date:**

Tue Oct 06 02:59:21 PM

**Date:**

4 Nov 2009

**Document's Title:**

fid

**Spectrum Title:**

User/Group Lars Troendlin/Pfaltz

LT445 nach Auf.

\_Proton CDCI3 uicon troend 14

**Frequency (MHz):**

(f1) 400.132

**Original Points Count:**

(f1) 24576

**Actual Points Count:**

(f1) 65536

**Acquisition Time (sec):**

(f1) 2.9688

**Spectral Width (ppm):**

(f1) 20.689

**Pulse Program:**

ZG30

**Temperature:**

300

**Number of Scans:**

16

**Acq. Date:**

Tue Oct 06 02:59:21 PM

