



Filename = 645 C-3.jdf
 Author = delta
 Experiment = single_pulse_dec
 Sample_id = S#412825
 Solvent = DMSO-D6
 Creation_time = 11-SEP-2009 11:09:42
 Revision_time = 25-SEP-2009 23:33:46
 Current_time = 25-SEP-2009 23:34:32

Comment = single pulse decouple
 Data_format = 1D_COMPLEX
 Dim_size = 26214
 Dim_title = 13C
 Dim_units = [ppm]
 Dimensions = X
 Site = ECX 400P
 Spectrometer = DELTA2_NMR

Field_strength = 9.2982153[T] (400[MHz])
 X_acq_duration = 1.048576[s]
 X_domain = 13C
 X_freq = 99.54517646[MHz]
 X_offset = 100[ppm]
 X_points = 32768
 X_prescans = 4
 X_resolution = 0.95367432[Hz]
 X_sweep = 31.25[kHz]
 Irr_domain = 1H
 Irr_freq = 395.88430144[MHz]
 Irr_offset = 5[ppm]
 Clipped = FALSE
 Mod_return = 1
 Scans = 140
 Total_scans = 140

X_90_width = 9.5[us]
 X_acq_time = 1.048576[s]
 X_angle = 30[deg]
 X_atn = 4.2[dB]
 X_pulse = 3.16666667[us]
 Irr_atn_dec = 21[dB]
 Irr_atn_noe = 21[dB]
 Irr_noise = WALTZ
 Decoupling = TRUE
 Initial_wait = 1[s]
 Noe = TRUE
 Noe_time = 2[s]
 Recvr_gain = 52
 Relaxation_delay = 2[s]
 Repetition_time = 3.048576[s]
 Temp_get = 22.3[dC]

